Instruction for Use EndoRotor® System

CE
2797

ESS-CONSOLE
PED-PURGE
EPR-TRAP
ESS-ROLLSTAND
ESS-VAC-DFLEX
3.2-PED-EGD-OP
3.2-PED-COL-OP
3.2-PED-EGD-F
3.2-EPR-EGD-OP
3.2-EPR-COL-OP
3.2-EPR-EGD-F
3.2-EPR-COL-O
EPR-FILTER
6.0-PED-EGD
6.0-EPR-EGD

EndoRotor® EndoSurgery System - Console
EndoRotor® Purge Kit
EndoRotor® Specimen Trap
EndoRotor® EndoSurgery System - Rollstand
EndoRotor® EndoSurgery System - Flex Vacuum Pump
EndoRotor® 3.2 mm Powered Endoscopic Debridement (PED) Catheter – Gastroscope Olympus/Pentax
EndoRotor® 3.2 mm Powered Endoscopic Debridement (PED) Catheter – Colonoscope Long Olympus/Pentax
EndoRotor® 3.2 mm Powered Endoscopic Debridement (PED) Catheter – Gastroscope Fuji
EndoRotor® 3.2 mm Endoscopic Powered Resection (EPR) Catheter – Gastroscope Olympus/Pentax
EndoRotor® 3.2 mm Endoscopic Powered Resection (EPR) Catheter – Colonoscope Long Olympus/Pentax
EndoRotor® 3.2 mm Endoscopic Powered Resection (EPR) Catheter – Gastroscope Fuji
EndoRotor® 3.2 mm Endoscopic Powered Resection (EPR) Catheter – Colonoscope Short Olympus
EndoRotor® Filter Set
EndoRotor® 6.0 mm Powered Endoscopic Debridement (PED) Catheter
EndoRotor® 6.0 mm Endoscopic Powered Resection (EPR) Catheter

1 Intended Use

The EndoRotor System is indicated for alimentary tract resection/excision and direct endoscopic necrosectomy (DEN) procedures to treat walled-off necrosis (WON).

2 Device Description

The EndoRotor® is a powered debridement tool for the removal of alimentary tract mucosa. The system consists of a power console, foot control, specimen trap with pre-loaded filter and the applied part, which is the single-use catheter.

3 Instruction Manual

This instruction manual contains essential information on using this instrument safely and effectively. This instruction manual is provided with the purchase of the EndoRotor® Console. Before use, thoroughly review this manual and the manuals of all equipment, which will be used during the procedure and use the instrument as instructed. For further information regarding the use of this product or to report any problems, please contact Interscope® using the appropriate information provided on the product label contact information packaged with each device; or contact your local distributor.

The information contained in this document is accurate at time of publication. Interscope® reserves the right to make changes to the product described in this manual without notice and without incorporating those changes to products already sold.

4 User Qualifications

The operator of this instrument must be a physician or medical personnel under the supervision of a physician and must have received sufficient training in clinical endoscopic technique. This manual, therefore, does not explain or discuss clinical endoscopic procedures.

5 Definitions

WARNING

Indicates a potentially hazardous situation, which if not avoided could result in death or serious injury.

PRECAUTION

Indicates a potentially hazardous situation, which if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices or potential equipment damage

NOTE

Indicates additional helpful information.

6 Warnings and Precautions

Follow the warnings and precautions described below when handling this system.

CONTRAINDICATIONS

The EndoRotor® should not be used as a primary means for resection in naïve lesions that both display non-lifting sign and are suspected as cancerous.

System/Disposable Warnings	
W1	It is important that the EndoRotor® system operator be familiar with the system Instructions for Use its precautions, procedures and safety issues.
W2	Failure to properly secure the console to the EndoRotor Roll Stand could result in the falling due to instability.
W3	Do not store the sterile packages containing the catheter in a place where they will become damaged, wet or improperly sealed.
W4	Do not use an instrument after the expiration date displayed on the sterile package.
W5	Inspect the sterile package for tears, inadequate sealing or water damage. If the sterile package shows any irregularities, the sterile condition of the instrument may have been compromised. Do not use the instrument.
W6	The catheter is a single use device. Do not attempt to reuse and sterilize the catheter. Failure to do so could result in pathology risk, performance issues and safety risk to the patient.
W7	Damage or irregularities such as breaches or loose connections may compromise patient or operator safety, pose an infection control risk or cause patient injury.
W8	Failure to use personal protective equipment could lead to blood, skin irritation, mucous and other potentially infectious material being transferred from the patient.
W9	Failure to load the vacuum tubing correctly into the vacuum control valve prevents the user from controlling vacuum and could result in injury to the patient.
W10	Do not insert the catheter into the working channel while operational.
W11	Do not articulate the endoscope if you are unaware of the location of the EndoRotor cutting surface. This could cause patient injury such as punctures, hemorrhage or injury to non-target tissue from mechanical resection.
W12	Do not probe using the catheter while it is operational. This could cause patient injury such as punctures, hemorrhage or injury to non-target tissue from mechanical resection.
W13	During the procedure, if you determine that the instrument is not operating properly – due to breakage, failure, or disconnection from the vacuum control valve – immediately stop using the device and carefully remove it to avoid patient injury.
W14	Do not withdraw the catheter into the endoscope quickly. This could damage the instrument.
W15	After using the instrument, disconnect all connections and dispose of the catheter in an appropriate manner. Failure to do so could pose an infection control risk.
W16	Remove the filter from the trap and place in a specimen jar for histopathology. Not doing so could result in loss of specimen and pathology error.
W17	Disconnect power to the EndoRotor® console before cleaning to avoid electrical shock.
W18	To avoid the risk of electrical shock, achieve electrical grounding reliability with proper connections. Connect the EndoRotor® system to hospital grade receptacles only.


W19	If this equipment is operated in the presence of high levels of electromagnetic interference (EMI) or highly sensitive equipment, interference may be encountered and the user should take whatever steps are necessary to eliminate or reduce the source of the interference. Diminished performance may lengthen procedure time.
W20	Medical Electrical Equipment needs special precautions regarding Electromagnetic Compatibility (EMC) and needs to be installed and put into service according to the EMC information provided in this Guide.
W21	Portable and mobile RF communications equipment can affect Medical Electrical Equipment.
W22	No modification of the EndoRotor® Console is allowed. Qualified Interscope personnel only must perform all service. Only cables and accessories supplied with the EndoRotor® Console must be used. Other cables and accessories could negatively affect EMC performance.
W23	Repair and/or modification to the EndoRotor® system by anyone other than qualified service personnel may significantly compromise the console's ability to perform effectively and void the equipment warranty.
W24	Only use Interscope components with the system.
W25	Do not use this instrument for any purpose other than its intended use.
W26	Use EndoRotor® supplied vacuum. Failure to do so can result in using an undesired negative pressure, which can result in degraded performance or unintentional cutting of tissue.
W27	When using this instrument, mild to moderate bleeding can occur with greater risk for bleeding in larger resections. Use physician discretion to manage bleeding using standard techniques.
W28	Do not attempt to use any substitute accessories such as foot pedals. Failure to do so will result in non-compliance for electrical safety and could pose a risk to the user and equipment.
W29	Follow appropriate hospital policy for disposal of electrical equipment in compliance with the WEEE Directive or where possible return to original end manufacturer for disposal.
W30	To avoid the risk of electric shock, this equipment must only be connected to a supply main with protective earth.
W31	A wait time of at least 90 seconds should be observed before pressing the stand-by button after powering down the console before powering the system on again.

System/Disposable Precautions	
P1	Use the EndoRotor® Catheter and Console only in combination with one another.
P2	Do not attempt to use the catheter in working channels smaller than 3.2 mm. Doing so may result in damage to the endoscope.
P3	Maintain sufficient distance from the console so that the flexible outer tubing is as straight as possible with no kinks or coils. Failure to do so may affect performance. If performance is affected, then straighten the outer tubing and operate the device.
P4	Never use excessive force to introduce the catheter when resistance is encountered. This may damage the instrument or endoscope.
P5	Do not run the system without first inserting into endoscope and priming with irrigation, as this may affect performance.
P6	Do not attempt to operate the catheter beyond 60 degrees of endoscope articulation as this may damage the endoscope or catheter.
P7	Do not make any connection inspections when the catheter is operational. This may result in operator or assistant injury, such as punctured gloves, tissue cuts or tissue abrasions.
P8	Do not force the connection of the proximal housing into the female connection on the console. If the hub does not slide in smoothly, verify the lock is in the unlocked  position with the lever completely to the left. Failure to do so may result in damage to the console and/or catheter.
P9	Do not soak the foot pedal or console. Doing so may result in equipment damage.
P10	Do not stack the EndoRotor Console on other electrical equipment. The console was validated using an IV pole mount. Failure to do so may result in electrical interference with the equipment it is stacked upon.
P11	In Endoluminal procedures, powered rotation of the cutter greater than 12 minutes may impact the performance of the catheter. If performance is affected after 12 minutes of powered rotation of the cutter then a new catheter may be needed.
P12	In DEN procedures, powered rotation of the cutter greater than 30 minutes may impact the performance of the catheter. If performance is affected after 30 minutes of powered rotation of the cutter then a new catheter may be needed.
P13	Do not operate the catheter beyond the point indicated within this IFU. Failure to do so may result in patient injury.
P14	Failure to center a replacement filter may result in loss of a vacuum seal when reattaching the Specimen Trap resulting in poor equipment performance.
P15	When removing the catheter from the console during procedure end follow the removal steps in the proper sequence. Failure to do so may result in console contamination.
P16	Failure to properly load the catheter tubing into the pump can cause crimping and limit flow.

P17	Do not twist or turn the Foot Control plug while inserting or removing the plug from the Console receptacle or the plug and/or Console receptacle could be damaged
P18	When running the Catheter without irrigation using the irrigation flow switch do not operate the Catheter greater than 3 minutes without re-engaging irrigation for 1 minute or flushing the Catheter using the prime button. Failure to do so may result in degraded performance.
P19	Do not retract the catheter into the working channel of the endoscope while the catheter is operational. Doing so may result in degraded performance.
P20	During necrosectomy procedures insert the catheter while the endoscope is in the stomach or duodenum. To avoid potential damage to a LAMS retract the catheter into the working channel so the instrument is not exposed to avoid damage to any cleared endoscopic stent that allows passage of the endoscopic to facilitate debridement.
P21	LAMS are often occluded with necrotic content. The user may debride necrosis within the LAMS with limited torqueing of the endoscope. Users should avoid probing the stent in the perpendicular plane to avoid damage to the stent.
P22	LAMS dislodgement may occur during necrosectomy. Users should minimize the number of times the endoscope is completely removed from the collection and patient in order to reduce this risk.

7 Nomenclature


7.1 EndoRotor® Catheter

1. Rotation handle
2. Aspiration tube
3. Saline bag connection
4. Cutting tip
5. Proximal housing/console interface
6. Flexible outer tubing

Figure 1 EndoRotor® Catheter


7.2 EndoRotor® Console

1. Specimen trap holder
2. Indicator light
3. Speed control switch
4. Standby button
5. Vacuum control release button
6. Prime button
7. Irrigation on/off button
8. Catheter interface locking lever
9. Catheter interface
10. Vacuum control valve
11. Irrigation pump
12. Foot control interface

Figure 2 EndoRotor® Console

7.3 EndoRotor® Foot Control

1. Console interface
2. Cutting control pedal
3. Vacuum control pedal

Figure 3 EndoRotor® Foot Control

7.4 EndoRotor® Roll Stand

1. Locking wheels
2. Handle
3. Console mounting post
4. Saline bag hooks
5. Basket for Control Pedal storage

Figure 4 EndoRotor® Roll Stand

8 Console Installation

1. User qualifications
 - a. Interscope® recommends that only qualified Interscope® personnel or hospital personnel who have been trained by Interscope® perform installation of the EndoRotor® Console.
2. The rear panel has a potential equalization connector as an alternative to equalizing grounding, which is normally achieved with the proper power cord and mains plug.
 - a. This can also be used for facility clinical engineers to inspect the Console for electrical leakage testing per their facility requirements.
3. Interscope® recommends using the EndoRotor® Roll Stand to mount the EndoRotor® Console as it provides the necessary stability. Refer to the EndoRotor® Roll Stand Installation Guide for instructions on assembling the Roll Stand.
4. When mounting the Console onto the Roll Stand ensure the top of the Console does not exceed a maximum height of 56 inches (144 cm).
5. Ensure the available power supply conforms to the system specification.
6. Do not position the Console such that one is unable to easily disconnect mains power via the appliance coupler.
7. Plug the Console into the appropriate power supply.

WARNING

Only use the power cord supplied with the Console. Other power cords could negatively affect EMC performance.

8. Ensure power is operational by pressing the on button, which is on the rear panel (Figure 5).
9. Depress the Standby button (Figure 6).
10. Ensure the indicator light located below the Standby button is amber (Figure 6).
11. Connect a Catheter to the Console, prime the Catheter by pressing the prime button (Figure 2) and ensure the indicator light located below the Standby button is green. Refer to figure #6.

Figure 5 On button

Figure 6 Stand by Button and Indicator Light

12. Determine and set the vacuum pressure and flow rate:

Vacuum Start Ranges*			
			Flow Rate 50 or 60 lpm
EndoRotor Cutter Speed	Mucosectomy		Necrosectomy
	Non-scarred	Scarred	550 mmHg
	High	50-100 mmHg	100-225 mmHg
Low	50-100 mmHg	100-225 mmHg	
*For general reference only. User must verify preferred vacuum and other system settings.			

9 Specimen Trap/Filter Preparation

Specimen Trap and Filter

1. The EndoRotor® Specimen Traps are packaged separately from the Catheter.
2. The Specimen Traps are preloaded with a nylon micron Filter. Multiple Specimen Traps may be used during a procedure in order to segregate tissue from multiple sites.
 - a. The EndoRotor® Filter Set (200 per pack) is available for order, EPR-FILTER (Figure 7).
3. The EndoRotor® Specimen Trap is non-sterile. It must be placed into the Specimen Trap holder (Figure 2) and connected to the vacuum source prior to the start of the procedure.

Figure 7 Filter replacements

Exchanging Filters

1. During a procedure the operator may elect to exchange a Filter to separate samples from a second location.
2. Obtain a replacement Filter from the EndoRotor® Filter Set (Figure 7).
3. Ensure any specimen still in the Catheter is flushed into the trap using the prime button.
4. Remove the vacuum tubing from the lower side of the trap.
5. Hold the lower half of EndoRotor® Specimen Trap. Using a second hand, rotate the upper half counter clockwise until they are unlocked (Figure 8).
6. Separate the halves then remove the used Filter using any gauge needle to lift the Filter and place the entire Filter into the pathology container (Figure 9.)

7. Replace with the new Filter and ensure it is centered before resealing the Specimen Trap (Figure 10).

Figure 8 Counter rotate trap halves

Figure 9 Separate Trap halves

Figure 10 Obtain new filter

10 Catheter Selection Guide

The following table will aid the user in selecting the appropriate Catheter based on the type of procedure and Endoscope being used.

SKU #	Minimum Working Channel Diameter	Cutter	Description
3.2-EPR-COL-OP	3.2 mm	 EPR Tip	EndoRotor® 3.2 mm Endoscopic Powered Resection (EPR) Catheter – Colonoscope Long Olympus/Pentax
3.2-EPR-COL-O	3.2 mm		EndoRotor® 3.2 mm Endoscopic Powered Resection (EPR) Catheter – Colonoscope Short Olympus
3.2-EPR-EGD-OP	3.2 mm		EndoRotor® 3.2 mm Endoscopic Powered Resection (EPR) Catheter – Gastroscope Olympus/Pentax
3.2-EPR-EGD-F	3.2 mm		EndoRotor® 3.2 mm Endoscopic Powered Resection (EPR) Catheter – Gastroscope Fuji
6.0-EPR-EGD	6.0 mm		EndoRotor® 6.0 mm Endoscopic Powered Resection (EPR) Catheter
3.2-PED-COL-OP	3.2 mm	 PED Tip	EndoRotor® 3.2 mm Powered Endoscopic Debridement (PED) Catheter – Colonoscope Long Olympus/Pentax
3.2-PED-EGD-OP	3.2 mm		EndoRotor® 3.2 mm Powered Endoscopic Debridement (PED) Catheter – Gastroscope Olympus/Pentax
3.2-PED-EGD-F	3.2 mm		EndoRotor® 3.2 mm Powered Endoscopic Debridement (PED) Catheter – Gastroscope Fuji
6.0-PED-EGD	6.0 mm		EndoRotor® 6.0 mm Powered Endoscopic Debridement (PED) Catheter

11 Catheter Inspection

1. When inspecting the instrument always wear appropriate personal protective equipment such as eye wear, a face mask, moisture resistant clothing and chemical resistant gloves that fit properly and are long enough so that skin is not exposed.
2. Before use, inspect the entire length of the Catheter for tears, cuts, peeling or other damage and do not use the instrument if any irregularity is found.
3. Ensure there are no breaches or holes in either the vacuum tubing or irrigation tubing.
4. While placing the insertion portion into the endoscope ensure there are no kinks in the tubing and run your fingertips over the entire length of the distal braided length of the Catheter during insertion to check for any crushed areas, broken areas or other damages. Do not use if such damage is found.

12 Catheter and Console Preparation for Use

1. Place the Foot Control in the operator's preferred position.
2. Connect the plug for the Foot Control to the lower left corresponding Console connection as shown in Figure 12. Ensure that the red dots on the Console and the plug are aligned prior to inserting the plug (Figure 13).

Figure 12 Foot Control Connection

Figure 13 Alignment of Red Dots

Figure 14 Placement of Specimen Trap

3. Ensure the EndoRotor® Specimen Trap is in the holder (Figure 14) and rotated counter clockwise to lock it in place.
4. Ensure the Catheter interface locking lever is in the unlocked position (Figure 15).
5. Insert the Catheter's proximal housing into the Catheter interface on the Console and ensure the connection aligns. If required, rotate the Catheter proximal housing slightly until it engages. Continue to push the Catheter's proximal housing in while rotating the Catheter interface locking lever to the locked position (Figure 16.)

Figure 15 Inserting Catheter proximal housing

Figure 16 Locking Catheter interface locking lever

Figure 17 Specimen Trap connection tube

6. Gently pull on the Catheter proximal housing to ensure it is secure.
7. If the user is unable to fully engage the lock to the right, then repeat steps 4 and 5 to ensure the connection is properly aligned.
8. On the bottom surface of the Catheter's proximal housing are 2 connected tubing sets (Figure 17). The shorter set without the spike fitting is the Specimen Trap connection tube.
9. To load the tubing, press the Vacuum control release button located directly to the right of the standby button (Figure 18). This will open the valve for ~25 seconds to ease insertion of the tubing.
10. Use 2 hands to stretch a 2-4 inch segment of tubing (Figure 19).
11. Gently push the stretched segment into the valve until it is completely seated (Figure 20).
12. Hold the tubing in place until the valve closes. One should also be able to see the tubing is completely within the valve. If a portion of the tube is not in the valve (Figure 21) repeat the steps to correct it.

Figure 18 Vacuum control release button

Figure 19 Tubing Stretch

Figure 20 Tubing insertion

Figure 21 Correct and incorrect tube placement

Figure 22 Connecting Vacuum Tube

Figure 23 Connecting Vacuum Tube

13. Connect the vacuum tubing directly to the trap with a press fit by pushing the tubing onto the trap until it is firmly seated and does not easily pull off (Figure 22).
14. Connect procedure room vacuum to the waste or overflow container.
15. Connect the overflow container to the Specimen Trap (Figure 23).
16. Use the vented spike on the irrigation tubing to attach to a saline IV bag.
17. Open the irrigation pump using one hand to lift the irrigation pump hood until the pump is open (Figure 24).

Figure 24 Inserting Tubing

Figure 25 Opening Pump

Figure 26 Closing Pump

18. Place the irrigation tubing set on top of the rollers in the pump with the flow from left to right (Figure 25.)
19. Close the irrigation pump by using one hand to lower the hood down until the tubing is engaged (Figure 26).

PRECAUTION Failure to properly load the tubing can cause crimping and limit flow. (Figure 27 and 28).

Figure 27 Incorrect

Figure 28 Correct

20. Press the prime button . If bubbles in the saline bag are visible, then the irrigation tubing is loaded backwards. Remove the irrigation tubing and reinstall.
21. Confirm the system is primed by observing fluid at the tip of the Catheter. If fluid is not present the user can press the prime button as many times as required.
22. Once the system has been primed and all user connections are properly engaged the amber indicator light turns green and the system is ready for use.
23. Select and switch the cutter speed to either High or Low based on user preference. During learning of technique, the physician can determine the preferred setting by trying the different speed options.

Inspection of Console operation

1. Depress the orange Foot Control to confirm the vacuum control valve is working properly. With a depressed vacuum Foot Control, the operator should see fluid being collected in the waste container and in the vacuum tubing line.
2. Ensure neither pedal is depressed so that the Catheter is not in operation mode and insert the Catheter until the cutting surface is viewable on the monitor.

Inspection of Catheter operation

1. Depress the blue Foot Control to engage the Catheter rotation functionality.
2. Ensure that the inner Catheter is rotating within the outer Catheter.

13 Purge Kit

Interscope markets an EndoRotor® accessory that is used in place of the Specimen Trap during procedures where tissue collection is not required.

Refer the Purge Kit instructions for use for details pertaining to this accessory.

14 Use

Inserting the Catheter into the endoscope & orienting

1. Insert the Catheter slowly.
2. The Catheter's maximum protrusion from the distal face of the endoscope must be minimized. Users shall ensure that the black markings (figure 29) on the distal end of the Catheter are kept in clear view at all times in order to ensure that the Catheter protrusion is maintained within acceptable limits.

WARNING Catheter protrusion beyond the distal marking may result in cutting unintended tissue.

NOTE Depending on the path the endoscope takes through the anatomy, the user may encounter resistance while inserting the catheter in the areas where the scope insertion tube is not fully straight. This resistance is normal and to be expected. Continue to insert the Catheter slowly through the resistance until it is fully deployed.

3. Degraded Catheter performance may occur if the endoscope is articulated to a point where the internal working channel is impinged. If degraded performance is noted (by observing reduced rotation speed or failure to cut) then reduce the articulation on the scope to confirm function and troubleshoot.
4. The Catheter Rotation handle is adjacent to the working channel to allow for repositioning of the outer cutter (Figure 30).
5. Before resecting tissue, reposition the Rotation handle counter clockwise or clockwise to get to the desired cutting orientation.
6. When using the Rotation handle small movements between 10-30 degrees are recommended. The greater the number of bends in the scope the finer the rotations should be. The operator can make larger movements but should build inertia in the rotation and prior to the cutter completing the movement the operator should rotate slightly either counter clockwise or clockwise to gain greater position precision.


Figure 29 Catheter distal markings to aid user in maintaining acceptable catheter protrusion

Figure 30 Rotation handle while operating scope

Resection

1. Do not depress either Foot Control until you are ready to begin resection.
2. Position the Catheter using the Catheter Rotation handle to the desired position.
3. Always use the demarcated line on the cutting surface to know exactly where the cutting surface is at all times. This line indicates the cutter position is 180 degrees of the solid black line. There are also 2 hashed black lines which indicate positions 90 degrees either side of the cutting surface to allow for resection at a more acute angle. A circumferential line at the tip indicates the optimal cutting area of the cutter (Figure 31).
4. Do not reposition the Catheter while operational and cutting. This can apply additional torque to the Catheter which can degrade performance or cause a device failure.
5. Always position the Catheter directly above the target to be resected prior to cutting.
6. Engage the blue Foot Control to initiate cutter rotation.
7. Engage the orange Foot Control within approximately 15 seconds or operation will cease.
8. Use the scope to direct the Catheter with either applying torque to the endoscope to roll the cutter from left to right or vice versa. Alternatively slide the endoscope from proximal to distal cutting. Use the greater articulation knob of the endoscope to lift away from the target tissue after the Catheter removes a resection site. Tap and retracting away from the target allow for constant reassessment of the removed site.
9. Allow vacuum to facilitate tissue into the cutting window. The active cutter should not be allowed to rest on the tissue when the user does not intend to resect. Do not bluntly probe or push into an organ wall with a fully operational Catheter.
10. Use the vacuum Foot Control to suction any fluid from the field of view without the rotation pedal. Resections normally require intermittent runtime alternating between resecting tissue and the user observing the results.


Figure 31 Positioning lines

- PRECAUTION** In Endoluminal procedures, powered rotation of the cutter greater than 12 minutes may impact the performance of the catheter. If performance is affected after 12 minutes of powered rotation of the cutter then a new catheter may be needed. In DEN procedures, powered rotation of the cutter greater than 30 minutes may impact the performance of the catheter. If performance is affected after 30 minutes of powered rotation of the cutter then a new catheter may be needed.
11. Relative anatomy is subjected to patient positioning and gravity. If saline begins to pool on or near the resection site, Console irrigation may be disengaged using the irrigation button  on the Console.
 12. Users should periodically flush the Catheter using the prime button or turn irrigation back on at 3-minute intervals to ensure performance of the Catheter does not degrade.

Flushing the Catheter

1. The user may elect to flush the Catheter to ensure the specimen is collected or in the event of a blockage.
2. Press the vacuum control valve button  and prime button  simultaneously. The operator should repeat this step until clear fluid without blood or tissue is visible in the aspiration tubing.
3. The user may also use a syringe filled with air and attach it to the aspiration tubing by first removing the aspiration tubing from the Specimen Trap.
4. Next attach the syringe to the aspiration tubing and press the vacuum control valve button .
5. While the valve is open plunge the air into the tubing and observe any debris that exits the Catheter.
6. Completing this action is recommended while the Catheter is rotating to ensure the cutting window is not occluded. During this step the Catheter is not capable of cutting without the aid of aspiration.
7. Reattach the aspiration tubing to the Specimen Trap and continue with the resection.

Withdrawing the Catheter from the endoscope

1. Cease Catheter operation and safely remove the Catheter from the endoscope.
2. Remove the Catheter with one hand while using the other hand to hold the Catheter at the working channel port with a gauze pad to prevent splashing of fluids.
3. Use the far hand to remove the Catheter while keeping the gauze pad on the biopsy port.
4. Once the Catheter is removed, any fluid will be collected on the pad and not splash the user.
5. Press the vacuum control valve button .
6. Remove the aspiration tubing.
7. Open the irrigation pump and remove the irrigation tubing.
8. Unlock the Catheter housing  and remove the Catheter. When performing the sequence in this order the Catheter will force itself out of the Console.

15 Disposal

15.1 Procedure End and Catheter Disposal

1. Place the Catheter in an appropriate bin for removal.
2. Ensure to follow local facility procedures for disposal of human waste.

15.2 Shutdown Procedure

1. Press the front panel standby button to turn off the Console.

Note It is not necessary to turn off the rear “off”  at the power supply. The operator can press the power prior to detaching the power cord from the power source.

2. Unplug the Console to complete the mains disconnect via the mains plug.

15.3 Cleaning and Disinfecting Procedure - Console

1. Do not immerse the Console in any fluid.
2. Clean the Console and or the Foot Control using isopropyl alcohol or disinfecting wipes to remove any blood or debris after each procedure.
3. Ensure the Console and Foot Control is free of debris and dry after each procedure.
4. Store the Console for later use in a clean and dry environment and not in direct sunlight.

15.4 Console Disposal Procedure

1. If the Console requires disposal, please follow local institutional policies for safe disposal of electronic equipment or contact Interscope®.

16 Preventative Maintenance and Repairs

1. All repairs or preventative maintenance shall be performed by qualified Interscope® personnel.
2. Interscope recommends an annual inspection of the consoles basic functionality. This inspection may be performed by any personnel who have read the Instruction for Use and are familiar with the use of the console. The inspection process does not require tools or equipment and may be obtained by contacting Interscope customer service. may be obtained by contacting Interscope customer service.
3. In order to mitigate device failure the customer may optionally arrange for annual preventative maintenance directly with Interscope or through the authorized distributor.

17 Replacement Components

WARNING: Only cables and accessories supplied with the EndoRotor® Console must be used. Other cables and accessories could negatively affect EMC performance.

1. A replacement power cord may be obtained by contacting Interscope.
2. A replacement foot control may be obtained by contacting Interscope.

18 Console Technical Specifications

18.1 Standards

This manual contains specifications in compliance with the following standards. Only qualified Interscope® personnel are authorized to repair the EndoRotor® Console. If the user seeks additional information such as diagrams, components list, etc. please contact Interscope® directly.

Standard Number	Standard Title
CFR	US Code of Federal Regulations -
ANSI/AAMI ES60601-1: 2005+A1: 2012	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-1:2005+ A1: 2012	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-2-18:2009	Medical electrical equipment - Part 2-18: Particular requirements for the basic safety and essential performance of endoscopic equipment
EN 60601-1-2: 2007/ AC: 2010	Medical electrical equipment - Part 1-2: General Requirements for Safety - Collateral Standard: Electromagnetic Compatibility - Requirements and Tests
CSA C22.2 No.60601-1-08	Medical Electrical Equipment - Part 1: General Requirements for basic Safety and essential performance

18.2 Operating environment – The EndoRotor® System is designed for use in a procedure room or operating suite at temperatures from at -29°C to 60°C and 30% to 85% relative humidity.

18.3 Transport and Storage Environment –

1. Catheter: Store the sterile packaged Catheter in a clean and dry environment. Do not store the instrument in direct sunlight. Ensure surrounding objects during storage do not crush the instrument.
2. The EndoRotor® Console and Catheters are designed for storage and transport at -29°C to 60°C and 30% to 85% relative humidity.

Physical Dimensions

Size	D 10.125 in x W 7.5 in x H 13.25 in
Size (Metric)	D 25.72 cm x W 18.42 cm H 33.65 cm
Weight	13.6 pounds (6.17kg)

Electrical

Input Voltage	100V to 240V
Frequency	50 – 60 Hz
Power Consumption	.75 – 1.4 amps
Duty Cycle for applied part	Continuous

Fuses

Size: 5mm x 20mm	Voltage: 250V C
Current (US): 2A (2)	EU Fuse Current: 1A (2)
Response Time: Slow	Breaking capacity: 150A

Essential Performance

Note This device does not provide Essential performance. If the device were to fail it would result in a delay of therapy. Alternative techniques can be used in its stead.

- Users may verify the following parameters are active at any time prior to the procedure.

Power Function

System operating Speeds

High (1750 revolutions per minute)

Low (1000 revolutions per minute)

Pinch Valve Release

Active/Inactive

Irrigation Prime

Produce Maximum Flow Rate

Indicator Light

Concurs set up steps have been executed

Electrical Safety Tables

Part I

Guidance and manufacturer's declaration - electromagnetic emissions		
<p>The EndoRotor® is intended for use in the electromagnetic environment specified below. The customer or the user of the EndoRotor® should assure that it is used in such an environment.</p> <p>NOTE: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.</p>		
Emissions test	Compliance	Electromagnetic environment — guidance
RF emissions CISPR 11	Group 1	<p>The EndoRotor® uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.</p> <p>The EndoRotor® is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.</p>
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Part II

Guidance and manufacturer's declaration — electromagnetic immunity			
The EndoRotor® is intended for use in the electromagnetic environment specified below. The customer or the user of the EndoRotor® should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment, If the user of the EndoRotor® requires continued operation during power mains interruptions, it is recommended that the EndoRotor® be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Part III

Guidance and manufacturer's declaration — electromagnetic immunity

The EndoRotor® is intended for use in the electromagnetic environment specified below. The customer or the user of the EndoRotor® should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the EndoRotor®, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 80 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the EndoRotor is used exceeds the applicable RF compliance level above, the EndoRotor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the [ME EQUIPMENT or ME SYSTEM].

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Part IV

Recommended separation distances between portable and mobile RF communications equipment and the EndoRotor®

The EndoRotor® is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the EndoRotor® can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the EndoRotor® as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.


NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

19 Manufacturer Information

Interscope, Incorporated
200 Commerce Drive
Northbridge, MA 01534
USA
Tel.: +1-800-461-4289
0031-(0)-73-3032541
www.interscopemed.com

Emergo Europe
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Limited Warranty

1. This Limited Warranty provides the following assurance for the customer who purchases an EndoRotor® system. This Limited Warranty is extended only to the buyer purchasing the System directly from Interscope® or from its authorized distributor or representative. The EndoRotor® system includes the Console and foot control, as well as the EndoRotor® Catheter, EndoRotor® Specimen Trap and Filter Set (hereinafter referred to as Single Use Components) and jointly referred to as the EndoRotor® Mucosal Resection System, unless specifically noted.
 - i. Should a System Component fail to function to Interscope's® published specifications during the term of this Limited Warranty (one year from the date of sale of a new System) Interscope® will either repair or replace the Console Component or any portion thereof.
 - ii. Should a Single Use Component fail to function to Interscope's® published specifications prior to its "use by" date Interscope® will replace the Single Use Component.
2. To qualify for this Limited Warranty, the following conditions must be met:
 - i. The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
 - ii. The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
 - iii. Interscope® must be notified in writing within thirty (30) days following discovery of a defect.
 - iv. The Product must be returned to Interscope® within thirty (30) days of Interscope® receiving notice as provided for above.
 - v. Upon examination of the Product by Interscope®, Interscope® shall have determined that: (i) the Product was not repaired or altered by any- one other than Interscope or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed annual maintenance services, if applicable, have been performed on the Product.
3. This Limited Warranty is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Interscope® be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the EndoRotor® Mucosal Resection System, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.
4. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this Limited Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Limited Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid.

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Standardnummer	Standardtittel
CFR	US Code of Federal Regulations -
ANSI/AAMI ES60601-1: 2005+A1: 2012	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-1:2005+ A1: 2012	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-2-18:2009	Medical electrical equipment - Part 2-18: Particular requirements for the basic safety and essential performance of endoscopic equipment
EN 60601-1-2: 2007/ AC: 2010	Medical electrical equipment - Part 1-2: General Requirements for Safety - Collateral Standard: Electromagnetic Compatibility - Requirements and Tests
CSA C22.2 No.60601-1-08	Medical Electrical Equipment - Part 1: General Requirements for basic Safety and essential performance

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