Etac Supply Center AB	Document type Clinical evaluation	Page 1(7)
Document title Bathroom devices clinical evaluation	Document ID 1A13408	Revision 2

Division Homecare

Clinical evaluation

Bathroom devices

**Bathboard Fresh, bathboard Rufus Plus, shower stools
Edge, Smart and Easy, toilet raiser Cloo, toilet raiser Hi-Loo,
toilet armrest Supporter, wall mounted shower seat Relax**

Responsible Dept	Replaces Document	Date (Valid from)
QA	N/A	2015-10-29
Issued		
		2015-11-02 Lars-Erik Gårdeskog, Quality manager
Reviewed		
		2015-11-02 Christina Ripe, international product manager, Occupational Therapist
Approved		
		2015-11-02 Lars-Erik Gårdeskog, Quality manager

Revision	Description	Date	Sign
1	First release	2015-10-29	L-E G
2	Editorial changes, medical recommendations added	2015-11-02	L-E G

Etac Supply Center AB	Document type Clinical evaluation	Page 2(7)
Document title Bathroom devices clinical evaluation	Document ID 1A13408	Revision 2

Table of Contents

1. Purpose and scope	3
2. References and definitions.....	3
3. Product description and intended use	3
4. Risks related to clinical use.....	4
5. Relevant essential requirements of MDD	4
6. Justification of method for clinical review	4
7. Objectives of the review and search strategy	5
8. Search of adverse events.....	6
9. Related Etac products	6
10. Final conclusion	6

Etac Supply Center AB	Document type Clinical evaluation	Page 3(7)
Document title Bathroom devices clinical evaluation	Document ID 1A13408	Revision 2

1. Purpose and scope

Clinical evaluation of bathroom products:

- Fresh, bathboard
- Rufus Plus, bathboard
- Edge, shower stool
- Smart, shower stool
- Easy, shower stool
- Cloo, toilet raiser
- Hi-Loo, toilet raiser
- Supporter, toilet armrest
- Relax, wall mounted shower seat

2. References and definitions

- The council of the European communities, council directive 93/42/EEC, medical devices (hereinafter MDD)
- ISO 13485:2003 Medical devices – Quality management systems – Requirements for regulatory purposes.
- ISO 14971:2012 Application of risk management to medical devices
- EN 12182:2012 Assistive products for persons with disability – General requirements and test methods

3. Product description and intended use

Fresh: Plastic bath board, with non-slip pattern and sturdy brackets to make washing safe and easy. The extra wide seating surface on one side makes it easier to get into and out of the bath tub.

Rufus Plus: Etac Rufus Plus bath board has a classic design and a capacity of 200 kg (441 lbs). The embedded aluminum tubes ensure the stability and weight capacity. The edge stops are easy to adjust for different bath widths.

Edge: three-legged shower stool.

Smart: four-legged shower stool.

Easy: four legged shower stool with round seat.

Intended use: Facilitates sitting when carrying out intimate hygiene.

Cloo: Height-adjustable toilet seat raiser that is fixed securely to the toilet. It is supplied with a 10 cm elevation. When armrests are included in the kit, these are to be fitted when mounting the product to the toilet. Cloo can also be fitted at an angle.

Etac Supply Center AB	Document type Clinical evaluation	Page 4(7)
Document title Bathroom devices clinical evaluation	Document ID 1A13408	Revision 2

Hi-Loo: A range of ergonomic toilet seat raisers that facilitates sitting down and getting up from the toilet. Hi-Loo can be used with or without armrests. Armrests provide greater security when seated and support when you sit down/get up from the toilet. The armrests can be flipped up, but can also be fixed in the lowered position.

Supporter: A range of toilet armrests that facilitates sitting down and getting up from the toilet. The armrests provide greater security when seated and support when you sit down/get up from the toilet. The armrests can be flipped up, but can also be fixed in the lowered position.

Relax: Etac Relax is a smart and elegant shower seat. It is perfect when space is limited – as it folds flat to the wall. You can choose the version that best suits your needs as Relax comes in many configurations and as easy to fit in any shower area. Scandinavian minimalistic design gives you maximum comfort and a more stylish bathroom.

Medical recommendations:

The product may be used in the following cases:

- Limited overall mobility
- Weakness of the joints and leg muscles
- Muscle diseases
- Total endoprosthesis of the hip joint

4. Risks related to clinical use

Risk analysis is performed for each product, for its intended use as well as for the risk of contributing to injuries. Remaining risks are minimized and found acceptable.

5. Relevant essential requirements of MDD

Relevant requirements in MDD are Annex 1.6a.

6. Justification of method for clinical review

There are no known published clinical investigations of this kind of products. However, there are studies showing the need of aids for bathing/showering/toileting to prevent falls.

Ergonomics at Home - Design for Safe Living and Home Care Doctoral Thesis Jenny Hjalmarson School of Technology and Health Royal Institute of Technology 2014

“It is estimated that about one third of all people over 65, and half of all people in the group over 80 years, fall once a year [69-71]. Old people mainly fall in their own home on a flat floor; most of the time falls are in the living room, bedroom or the hallway [69, 72, 73]. Elderly people often experience a fear of falling [74-77]; the fear is also known to increase the risk of falling [78, 79]. Earlier studies show that it takes longer time to get up after a fall the older you get [80].

Etac Supply Center AB	Document type Clinical evaluation	Page 5(7)
Document title Bathroom devices clinical evaluation	Document ID 1A13408	Revision 2

The cost of injuries caused by people falling is high and it will increase when the group of elderly increases [81]. Injuries can be directly related to the specific fall, but even a seemingly minor fall can lead to fatal consequences for the person if he or she is not able to get up soon enough [82]. P 15-16

" Kroemer describes how to design for older people; his examples regarding the bathroom are simple but clear; "Provide ample space for dressing, washing, toileting, especially when less mobile because of age or disability; grab bars may be desirable" P 14.

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Authors: Inga-Lill Boman, Louise Nygård, Lena Borell

Sektionen för arbetsterapi, NVS, Karolinska Institutet, released by Hjälpmedelsinstitutet, have motivated the need for aids in the bathroom to support the user or shortening the distance when seating or raising.

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"Every year almost 75000 people need hospital stay after falling in the bathroom. Almost 50000 of them are 65 years or older."

"It is easy to slip on wet floors. Antislip devices for bathrooms are nowadays available in a great variety of colors and contours. A grab rail makes it easier to get out of the bath tub and can be an extra safety in the shower. Keep soap, schampoo and towel within reachable distance."

(Free translation from swedish original)

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In bathrooms it is most common with injuries from contact with toilet porcelain (38%), bath tubs (30%) and showers (15%). Fall accidents is the dominating cause and stands for 85% of injuries in cintact with toilet porcelain, 75% of accidents with bath tubs and 70% of accidents in showers.

(Free translation from swedish original)

7. Objectives of the review and search strategy

The objectives are to estimate the patient risk of the product by searching existing clinical studies performed on similar products and evaluate whether this data is sufficiently enough to sustain as a clinical evaluation without the need of a clinical study of the product.

Etac Supply Center AB	Document type Clinical evaluation	Page 6(7)
Document title Bathroom devices clinical evaluation	Document ID 1A13408	Revision 2

8. Search of adverse events

Incident report search

FDA MAUDE database:

Product class: Adaptor hygiene. Event type: Death. Date Report Received By FDA: 01/01/2013 – 30/09/2015

Result: Four reported deaths where medical devices were involved. Two cases were caused by collapsing seats, one user slipped off the seat, and one fell despite holding on to a grab rail.

Product class: Adaptor hygiene. Event type: Injury. Date Report Received By FDA: 01/01/2013 – 30/09/2015

Result: 59 reports related to:

- Collapsed seat
- Seat cracked
- Collapsing leg
- Wheel coming loose

MHRA, safety warnings, alerts and recalls:

Alert type: Medical device alert, and Field safety notice

Issued after 01/01/2013

Issued before 30/09/2015

Result: No information of safety alerts or recalls related to bathroom medical devices.

9. Related Etac products

- Etac has a wide range of similar products which have been in the assortment for more than 30 years.

We have with our experience of developing and manufacturing similar products a very good platform for that the bathroom products are well motivated and that they in a safe way compensates for impairments.

10. Final conclusion

We conclude that:

- There is no known clinical investigation published within the product range
- There are no known registered adverse events on similar products where the product function is questioned
- There are studies showing that manual medical devices in the bathroom helps to avoid fall accidents

Etac Supply Center AB	Document type Clinical evaluation	Page 7(7)
Document title Bathroom devices clinical evaluation	Document ID 1A13408	Revision 2

- The products have been evaluated through test methods required in the harmonized standard EN 12182
- Etac has more than 30 year experience of development, manufacturing and distribution of manual medical devices for bathrooms.

The conclusion is that the current information has addressed all the concerns enough. The product is deemed safe enough to be used for the previously explained intended use. As a result, no clinical investigation is needed. Conformance with the medical device directive has thus been demonstrated.

Annex 1:

Lars-Erik Gårdeskog has been responsible for product risk management process at Etac Supply Center AB since 1993.

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Lars-Erik Gårdeskog zodpovídal za proces řízení rizik výrobků v ETAC Supply Center AB od roku 1993.