EN Fermathron® plus

Sodium hyaluronate 30 mg/2,0 ml

Intended Use
For the relief of pain and stiffness of the knee, hip and ankle joints in patients with mild to moderate osteoarthritis resulting from degenerative and traumatic changes to the synovial joint.

Fermathron® plus

Presentation
Fermathron® plus is a clear solution of sterile 1.5% HMW sodium hyaluronate in phosphate buffered saline contained in a pre-filled syringe for single intra-articular injection into the synovial space of the joint.

2,0 ml of Fermathron® plus, sterilised by filtration, is enclosed within a glass, ready to use, disposable syringe. The syringe is packed within a blister pack and an outer cardboard carton. The outer surface of the pre-filled syringe has been sterilised by ethylene oxide.

Description and composition

Sodium hyaluronate is a long chain polysaccharide made up of repeating disaccharide units of D-glucuronic acid and N-acetyl-D-glucosamine, which occurs naturally in the human body. Fermathron® plus has a pH similar to synovial fluid.

Each ml of solution contains: Sodium hyaluronate 15 mg, sodium chloride Ph.Eur. 9,0 mg, Disodium hydrogen phosphate Ph.Eur. 0,27 mg, Sodium dihydrogen phosphate Ph.Eur. 0,092 mg, Water for injections Ph.Eur. q.s.

Dosage and Administration

Injection of Fermathron® plus should only be carried out by a Healthcare Professional trained in the technique.

The recommended dosage regimen for patients with mild to moderate osteoarthritis of the knee is up to three weekly injections of 2,0 ml into the synovial space of the knee joint. The dosage regimen should be adapted by the Healthcare Professional for injection into the synovial space of the hip and ankle joints. It is recommended that injections in the hip and ankle joints are performed using ultrasound or fluoroscopic guidance.

Clean the skin around the injection site with antiseptic and allow to dry before the injection is given.

If joint effusion is present it should be aspirated before injection of Fermathron® plus. The contents of the syringe are sterile and should be injected using a sterile needle of an appropriate size (18 to 21 gauge needle is recommended). The syringe is fitted with a Luer Lok™ (6%) to ensure secure attachment of the needle to the syringe.

Discard the syringe and needle after single use.

Uses

For the relief of pain and stiffness of the knee, hip and ankle joints in patients with mild to moderate osteoarthritis resulting from degenerative and traumatic changes to the synovial joint.

The duration of effect in patients with mild to moderate osteoarthritis of the knee joint is up to six months. Duration of effect in the hip and ankle joints has not been demonstrated.

The performance of Fermathron® plus is due to its biocompatibility and physicochemical properties. The HMW sodium hyaluronate contained in Fermathron® plus is a biopolymer composed of repeating disaccharide units of D-glucuronic acid and N-acetyl-D-glucosamine and though it is biosynthesised by the bacterium *Streptococcus equi* it has been shown to be the same as the sodium hyaluronate which is found in the human body. The hyaluronan supplements the hyaluronan found naturally in the synovium but which has been depleted by degenerative and traumatic changes to the synovial joint.

Contraindications
Do not inject Fermathron® plus if the area of the injection is infected or where there is evidence of skin disease. Do not use in patients with known hypersensitivity to sodium hyaluronate.

Special warnings and precautions
Do not use if the sterile packaging has been damaged.

Do not use after the expiry date. Sodium hyaluronate is manufactured by fermentation of *Streptococcus equi* and rigorously purified. However, the physician should consider that the immunological and other potential risks that can be associated with the injection of any biological material.

There is no evidence concerning the safety of Fermathron® plus in human pregnancy and lactation. Administration during pregnancy and lactation is at the discretion of the orthopaedic surgeon.

Follow national or local guidelines for the safe use and disposal of needles. Obtain prompt medical attention if injury occurs.

Fermathron® plus is a single use device and is intended to be used for a single patient only. If Fermathron® plus is used for a second patient the sterility will be compromised and there is a risk of foreign body reaction and/or infection. Fermathron® plus should not be re-sterilised as the device performance may be compromised. Fermathron® plus should be used with a sterile needle that should be discarded after single use.

Adverse reactions

Transient pain and swelling may occur following intra-articular injections.

Rarely an inflammatory reaction, septic arthritis, or arthralgia could occur which may or may not be associated with Fermathron® plus.

Incompatibilities

Fermathron® plus has not been tested for compatibility with other substances for intra-articular injection. Therefore the mixing or simultaneous administration with other intra-articular injectables is not recommended.

Storage

Store between 2°C and 25°C. Do not freeze. Protect from light. Sterile product for single use only. Do not use after expiry date.

Symbols Used on the Packaging	
LOT	Batch code.
REF	Catalogue number.
	Use-by date.
	Do not re-use.
	Temperature limit (store between 2°C and 25°C)
	Keep away from sunlight.
	Do not freeze.
	Do not use if the packaging is damaged.
STERILE A	Sterilised using aseptic processing techniques.
STERILE EO	Sterilised using ethylene oxide.
	Date of manufacture.
	Manufacturer.
	Distributor.
	Consult instructions for use.

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