# **PHASIX™ Mesh**

Fully Resorbable Implant  
for Soft Tissue Reconstruction

## Instructions for Use

Mode d'emploi • Gebrauchsanweisung  
Istruzioni per l'uso • Instrucciones de uso  
Gebruiksaanwijzing • Instruções de utilização  
Οδηγίες Χρήσης • Brugervejledning  
Bruksanvisning • Käyttöohjeet  
Bruksanvisning • Instrukcja użycia  
Használati útmutató • Návod k použití  
Kullanım Talimatları • Инструкции по применению  
使用説明 • 사용 설명서

**STERILE EO**

Single Use

Do Not Resterilize

Absorbable

Only

Do not use if package is damaged

Manufacturer:

Davol Inc.  
Subsidiary of C. R. Bard, Inc.  
100 Crossings Boulevard  
Warwick, RI 02886 USA  
1-401-825-8300 • 1-800-556-6275

**EC REP**

BD Switzerland Sàrl  
Terre Bonne Park – A4  
Route de Crassier 17  
1262 Eysins, Switzerland

**CE**  
2797

**Medical Services & Support**  
**Clinical Information Line**  
**1-800-562-0027**

PK3799455  
REV. 2020/03

**BAARD**

DAVOL INC.

## PHASIX™ Mesh

Fully Resorbable Implant  
for Soft Tissue Reconstruction

### PRODUCT DESCRIPTION

PHASIX™ Mesh is a fully resorbable synthetic mesh implant prepared from poly-4-hydroxybutyrate (P4HB). P4HB is produced from a naturally occurring monomer and is processed into monofilament fibers and then knitted into a surgical mesh. PHASIX™ Mesh degrades through a process of hydrolysis and a hydrolytic enzymatic digestive process. It has been developed to enforce areas where weakness exists while minimizing the variability of resorption rate (loss of mass) and strength to provide support throughout the expected period of healing.

Preclinical implantation studies indicate that PHASIX™ Mesh retains approximately 70% of its strength at 12 weeks. Absorption of the mesh material will be essentially complete within 12 to 18 months.

### INDICATIONS

PHASIX™ Mesh is indicated to reinforce soft tissue where weakness exists in patients undergoing abdominal, plastic, and reconstructive surgery in ventral hernia repair and other abdominal fascial defect procedures.

### CONTRAINDICATIONS

1. Because PHASIX™ Mesh is fully resorbable, it should not be used in repairs where permanent wound or organ support from the mesh is required.
2. PHASIX™ Mesh is contraindicated for use in the repair of pelvic organ prolapse.
3. PHASIX™ Mesh is contraindicated for use in the treatment of stress urinary incontinence.

### WARNINGS

1. PHASIX™ Mesh must not be put in direct contact with bowel or viscera.
2. Device manufacture involves exposure to tetracycline hydrochloride and kanamycin sulfate. The safety and product use for patients with hypersensitivities to these antibiotics is unknown. Use of this device in susceptible patients with known allergies to tetracycline hydrochloride or kanamycin sulfate should be avoided.
3. The safety and effectiveness of PHASIX™ Mesh in the following applications has not been evaluated or established:
  - a. Pregnant women
  - b. Pediatric use
  - c. Neural and cardiovascular tissue
4. If an infection develops, treat the infection aggressively. Consideration should be given regarding the need to remove the mesh. An unresolved infection may require removal of the device.
5. To prevent recurrences when repairing hernias, the PHASIX™ Mesh must be large enough to provide sufficient overlap beyond the margins of the defect. Careful attention to mesh fixation placement and spacing will help prevent excessive tension or gap formation between the mesh and fascial tissue.
6. Prior to use, carefully examine package and product to verify neither is damaged and that all seals are intact. Do not use if the foil pouch or package is damaged or open, or if the center of the temperature indicator on the foil pouch is black.
7. This mesh has been designed for single use only. Reuse, reprocessing, resterilization or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the mesh and may lead to mesh failure which may result in injury to the patient.  
Reuse, reprocessing, resterilization or repackaging may also create a risk of contamination of the mesh and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the mesh may lead to injury, illness or death of the patient or end user.

8. If unused PHASIX™ has been in contact with instruments or supplies used on a patient or contaminated with body fluids, handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

### PRECAUTIONS

1. Please read all instructions prior to use.
2. Only physicians qualified in the appropriate surgical techniques should use this device. Users should be familiar with mesh strength and size requirements. Improper selection, placement, positioning and fixation of the device can cause subsequent undesirable results.

### ADVERSE REACTIONS

In preclinical testing, PHASIX™ Mesh elicited a minimal tissue reaction characteristic of foreign body response to a substance. The tissue reaction resolved as the mesh was resorbed. Possible complications may include, but are not limited to infection, seroma, pain, mesh migration, wound dehiscence, hemorrhage, adhesions, hematoma, inflammation, allergic reaction, extrusion, erosion, fistula formation and recurrence of the hernia or soft tissue defect.

### DIRECTIONS FOR USE

PHASIX™ Mesh may be cut to the shape or size desired for each specific application. To prevent recurrences when repairing hernias, a mesh larger than the defect is required to ensure adequate coverage. The mesh is to be positioned so that its edges extend beyond the margins of the defect.

It is recommended that surgical fixation be performed using absorbable or permanent sutures, or absorbable or permanent tackers placed 6 to 12 mm (¼ to ½ inches) apart at a distance approximately 6 mm (¼ inch) from the edge of the mesh.

### HOW SUPPLIED

PHASIX™ Mesh is available in sterile packages as a single flat mesh of varying widths and lengths.

### STORAGE

Store at room temperature. Avoid prolonged exposure to elevated temperatures. If the center of the temperature indicator on the box is BLACK, check the temperature indicator on the foil pouch. If the center of the temperature indicator on the foil pouch is BLACK, DO NOT USE PRODUCT.

### TRACEABILITY

A traceability label that identifies the type, size, expiration date and lot number of the device is attached to every package. This label should be affixed to the patient's permanent medical record to clearly identify the device that was implanted.

If you experience a product failure, please contact Davol, Inc. at 1-800-556-6275 for instructions on returning the product.

Copyright © 2010 - 2020 C. R. Bard, Inc. All Rights Reserved. Bard, Davol and Phasix are trademarks and/or registered trademarks of C. R. Bard, Inc.

	U. S. Federal law restricts this device to sale by or on the order of a physician.		
	Do not use if package is damaged		
	Contents		Actual Size
	Absorbable		Square
	Circle		Rectangle
	Do not use if the center of the temperature indicator is black.		

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**Document Number:** BAW3791660

**Revision:** 1

**Title:** Phasix™ Mesh, CE-marked, m/I IFU Graphics

All dates and times are in Eastern Standard Time.

### DA CR-15965 TGatta Phasix IFUs

#### DC Change Request

Name/Signature	Title	Date	Meaning/Reason
Denise Menard (DMENARD)	SUPERVISOR, DOCUMENTATION CTRL		
Documentation Placeholder (DOCHOLDER)	Place Holder		
Julia Kue (JKUE)	DOCUMENTATION SPECIALIST		
Lester Torrey (LTORREY)	ASSOCIATE DIRECTOR, QA		
Michelle DAmbra (10300496)	Documentation Clerk	08 Apr 2020, 07:52:15 AM	Approved

#### Team Review/Approval

Name/Signature	Title	Date	Meaning/Reason
Tony Rennick (10308044)	Packaging Coordinator Temp	08 Apr 2020, 03:29:08 PM	Approved
Katherine Earle (KEARLE)	RA Reg Affairs Specialist	09 Apr 2020, 09:45:31 AM	Approved
Evans Kipyego (EKIPYEGO)	Sr. AQE	09 Apr 2020, 10:56:13 AM	Approved
Tara Gatta (TGATTA)	PKG PACKAGING COORDINATOR	09 Apr 2020, 11:23:44 AM	Approved

#### Release

Name/Signature	Title	Date	Meaning/Reason
Michelle DAmbra (10300496)	Documentation Clerk	09 Apr 2020, 12:21:44 PM	Approved

#### Notify

Name/Signature	Title	Date	Meaning/Reason
Kathy Strange (KSTRANGE)	RA Sr. Mgr. RA OPs/Int'l	09 Apr 2020, 12:21:45 PM	Email Sent
Juan DeCardenas (JDECARDENAS)	PD Director PD/AME	09 Apr 2020, 12:21:45 PM	Email Sent
Pam Fotopoulos (PFOTOPOULOS)	PD - Project Mgr Labeling	09 Apr 2020, 12:21:45 PM	Email Sent
Sonia Rodriguez (SRODRIGUEZ)	SR STRATEGIC SOURCING ANALYST	09 Apr 2020, 12:21:45 PM	Email Sent
Tara Gatta (TGATTA)	PKG PACKAGING COORDINATOR	09 Apr 2020, 12:21:45 PM	Email Sent

#### Quick Approval

**Approve Now**

Name/Signature	Title	Date	Meaning/Reason
Michelle DAmbra (10300496)	Documentation Clerk	22 Apr 2020, 06:59:15 AM	Approved

Released

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Copyright © 2010 - 2020 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2010 - 2020 C. R. Bard, Inc. Alle Rechte vorbehalten. Bard, Davol und Phasix sind Marken und/oder eingetragene Marken von C. R. Bard, Inc.
Copyright © 2010 - 2020 C. R. Bard, Inc. Tutti i diritti riservati. Bard, Davol e Phasix sono marchi commerciali e/o registrati di C. R. Bard, Inc.
Copyright © 2010 - 2020 C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2010 - 2020 C. R. Bard, Inc. Alle rechten voorbehouden.
Copyright © 2010 - 2020 C. R. Bard, Inc. Todos os direitos reservados.